


Angiography Drapes With Multiple Fenestration 2



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Description	Qty. / Ctn.
Cat. Code	Angiography drape	15pcs
AGD-103H		

Drape features:

- Absorbent impervious fabric (210x330cm)
- Single clear panel sheet
- Dual oval Apertures with adhesive (7x9cm)
- Single Aperture with adhesive (5x7cm)