


Cleft Palate Drape 3



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	CPD-002	Contents:	18pcs
		4 Tape Strips (5x15cm)	
		1 Headbar Drape (106x81cm)	
		1 Cleft Palate Drape (180x274cm)	

Drape features:

- Turban/head drape-Adhesive
- Cleft Palate-Fenestration (10x12cm) / Adhesive fenestration (1x6cm)