


Neuro Drape with IGS cover



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	NSD-002-IGS	Neuro Drape with IGS covers	10pcs

Drape features:

- Absorbent impervious fabric (46x72cm)
- Attached with tying string fluid collection pouch and IGS covers
- Fenestration with full incise (18x27cm)
- Cord/tube holders