


Shoulder Drape (Beach Chair Position)



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	SDD-001	Shoulder Drape (Beach Chair Position)	8pcs

Drape features:

- Attached with tying string fluid collection pouch
- Adhesive-backed fenestration
- Elastic fenestration (13.5x14.5cm)
- Flexible cord/tube holders